

Option 1

Attach above fittings directly to your current device tubing.

Option 2

Buy FlexiPort cuffs with tubing and connectors to fit your blood pressure devices without making changes.

No Tubes or Connectors Required

One-Tube

Two-Tube



Size (each)

Range

FlexiPort

Screw

Tri-Purpose

Locking

Bayonet

Screw

Locking

Tri-Purpose & Empty Tube

Inflation Sys.


Small Infant (6)	7-10 cm	REUSE-06	REUSE-06-1SC	REUSE-06-1TP	REUSE-06-1MQ	REUSE-06-1HP	n/a	n/a	n/a	n/a
Infant (7)	9-13 cm	REUSE-07	REUSE-07-1SC	REUSE-07-1TP	REUSE-07-1MQ	REUSE-07-1HP	REUSE-07-2SC	REUSE-07-2MQ	REUSE-07-2TP	REUSE-07-2BV
Small Child (8)	12-16 cm	REUSE-08	REUSE-08-1SC	REUSE-08-1TP	REUSE-08-1MQ	REUSE-08-1HP	REUSE-08-2SC	REUSE-08-2MQ	REUSE-08-2TP	REUSE-08-2BV
Child (9)	15-21 cm	REUSE-09	REUSE-09-1SC	REUSE-09-1TP	REUSE-09-1MQ	REUSE-09-1HP	REUSE-09-2SC	REUSE-09-2MQ	REUSE-09-2TP	REUSE-09-2BV
Small Adult (10)	20-26 cm	REUSE-10	REUSE-10-1SC	REUSE-10-1TP	REUSE-10-1MQ	REUSE-10-1HP	REUSE-10-2SC	REUSE-10-2MQ	REUSE-10-2TP	REUSE-10-2BV
Adult (11)	25-34 cm	REUSE-11	REUSE-11-1SC	REUSE-11-1TP	REUSE-11-1MQ	REUSE-11-1HP	REUSE-11-2SC	REUSE-11-2MQ	REUSE-11-2TP	REUSE-11-2BV
Adult Long (11L)	25-34 cm	REUSE-11L	REUSE-11L-1SC	REUSE-11L-1TP	REUSE-11L-1MQ	REUSE-11L-1HP	REUSE-11L-2SC	REUSE-11L-2MQ	REUSE-11L-2TP	REUSE-11L-2BV
Large Adult (12)	32-43 cm	REUSE-12	REUSE-12-1SC	REUSE-12-1TP	REUSE-12-1MQ	REUSE-12-1HP	REUSE-12-2SC	REUSE-12-2MQ	REUSE-12-2TP	REUSE-12-2BV
Large Adult Long (12L)	32-43 cm	REUSE-12L	REUSE-12L-1SC	REUSE-12L-1TP	REUSE-12L-1MQ	REUSE-12L-1HP	REUSE-12L-2SC	REUSE-12L-2MQ	REUSE-12L-2TP	REUSE-12L-2BV
Thigh (13)	40-55 cm	REUSE-13	REUSE-13-1SC	REUSE-13-1TP	REUSE-13-1MQ	REUSE-13-1HP	REUSE-13-2SC	REUSE-13-2MQ	REUSE-13-2TP	REUSE-13-2BV

FlexiPort™ Reusable Cuffs

True standardization in a snap.

Longlasting, easy to use, and affordable, Welch Allyn cuffs are the top choice for hospitals and clinics around the world. Why? Because Welch Allyn can offer you virtually every type, size, and configuration of cuffs, along with the ability to simplify inventory like no other cuff provider with our innovative FlexiPort™ connection, which makes every cuff a one- or two-tube cuff.

- Easy one-hand attachment/detachment of tubing from FlexiPort makes changing cuffs much faster and easier
- Tested to withstand demanding clinical environments
- 3-year warranty
- Meets all the latest clinical guidelines for proper fit from AAMI and AHA
- Standardizing can reduce part numbers by up to 60%




What's FlexiPort?

The FlexiPort fitting provides a single-point of connection between devices and cuffs, making any cuff a one- or two-tube cuff.



Patented FlexiPort™ technology allows one cuff to work in virtually any patient care area.



- Primary sizes are color-coded to make it easy to find the proper cuff size
Latex-free material minimizes risk of allergic reactions
- Rotatable port reduces stress to cuff tubing and port, improving patient comfort and cuff durability
- No sewn seams to tear or fray
Antimicrobial coating helps prevent the growth of mildew and bacteria
- Folded edge reduces risk of cuts and scrapes, maximizing patient comfort

The right cuff. The right tubes. At the right time.

The unique FlexiPort fitting enables true cuff standardization throughout your facility for the first time. With FlexiPort, finding the right cuff is made simpler than ever before by eliminating annoying hassles and frustrations of tube and connector configurations.

