

**COMBIMATT**

Double use vacuum/transport mattress



Finally the problem of stocking vacuum mattresses is solved. If used on one side the Combimatt is a soft mattress for stretchers, if used on the other side it is an excellent vacuum mattress with longitudinal canals containing the air chambers. Used as a mattress for stretcher, Combimatt offers characteristics that optimise the transportation of the patient, as it absorbs all mechanical vibrations and improves thermal isolation.

The on/off valve made of chromed brass with a built-in fitting for rubber tube is resistant to any temperature and guarantees unlimited duration.

Combimatt is supplied with

- three one-piece 50 mm nylon security belts with plastic buckle

Functioning temperature: from -10 to +50 °C

**TECHNICAL DATA**

Dimensions	2030 x 990 mm
Weight	6,4 kg
Covering superior part	Tarpaulin 600 D
Covering inferior part	Tarpaulin 700 D
Comfort system	High density closed cell expanded PE
Separate chambers	Yes
Internal material	Polystyrene
Handles	8
Aluminium one way valve	Yes
Space for spinal board	Yes
Distal tapering	Yes
Radio compatible	Yes
HFW	Yes
Compatible with traction systems	Yes
Colour	Orange / blue

Combimatt

QM22800A



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>